

The logo for Aquatic Diagnostics Ltd features the word "AQUATIC" in a bold, white, sans-serif font, set against a dark blue rectangular background. Below this, the words "Diagnostics Ltd" are written in a smaller, blue, sans-serif font. A light blue wavy graphic element is positioned behind the text.

**AQUATIC**  
Diagnostics Ltd

Aquatic Diagnostics Ltd.,  
Institute of Aquaculture,  
University of Stirling,  
Stirling, Scotland,  
FK9 4LA

Telephone: +44 (0) 1786 466568  
Fax: +44 (0) 1786 4672133  
E-mail: [aquaticdiagnostics@stir.ac.uk](mailto:aquaticdiagnostics@stir.ac.uk)  
<http://www.aquaticdiagnostics.com>

The logo for Aquatic Diagnostics Ltd features the word "AQUATIC" in a bold, white, sans-serif font, set against a dark blue rectangular background. Below this, the words "Diagnostics Ltd" are written in a smaller, blue, sans-serif font. A light blue wavy graphic element is positioned behind the text.

**AQUATIC**  
Diagnostics Ltd

*Anti-Infectious Salmon  
Anaemia Virus (ISAV)*  
monoclonal antibody

Product no: P10

The text "Product Information" is written in a large, white, sans-serif font. It is positioned over a light blue wavy graphic element that curves across the bottom of the page.

Product Information

## Product Description

The monoclonal antibody (MAb) against *Infectious Salmon Anaemia Virus (ISAV)* is specific for this virus. The specificity of the MAb has been tested against a range of viral pathogens that infect fish, including Noda virus and (Infectious pancreatic necrosis virus) IPN virus. The MAb is of an IgG1 isotype.

## Use of product

The MAb is recommended for use in immunohistochemistry (IHC), but can also be used in IFAT. The optimal conditions for use of this product vary depending on the procedure used. The user must determine the suitability of the product for a particular procedure. This product is for *in vitro* use only.

## Vial Contents

Each vial contains 200 µg of lyophilised protein prepared from bovine-free culture medium and contains no animal-derived stabilisers. This is sufficient for between 100-200 tests depending on the area of tissue to be screened in IHC.

The product should be reconstituted as follows:

- Add 1 ml of phosphate buffered saline (PBS) (see buffers) to the vial, then transfer the contents of the vial into 9 ml of PBS so that the total volume equals 10 ml.



## Certificate of Analysis

*Anti-Infectious Salmon Anaemia Virus (ISAV) monoclonal antibody*

Product no.

Batch no.

Date of expiry

### Activity in IHC

Cells infected with the virus appear golden brown in colour when stained with DAB.



## Storage

Store at  $-20^{\circ}\text{C}$  or below prior to reconstitution. For prolonged storage, the Mab solution should be stored at  $-20^{\circ}\text{C}$ , or below, as working aliquots. Repeated freeze thawing of the product should be avoided.

## Suggested protocol for the detection of ISAV in fixed tissue sections by immunohistochemistry

This procedure has been developed to work on tissues fixed in 10% buffered formalin for 24 hours. Individual protocols may have to be developed depending upon the tissue examined, fixation etc.

### Procedure (wear gloves)

- Prepare paraffin-embedded tissue sections.
- Dewax and rehydrate sections in xylene (2 x 5min), 100% ethanol (5 min), 70% ethanol (3 min), then rinse in distilled water.
- Place slides in a humid chamber.
- Keep sections moist at all times - do not allow them to dry out.
- Mark rings around the tissue sections using a wax PAP pen.
- Block endogenous peroxidase activity by incubating the slides for 10 min at room temperature ( $\approx 22^{\circ}\text{C}$ ) with  $\text{H}_2\text{O}_2$  in methanol (see buffers).
- Wash the slides three times with Tris buffered saline (TBS) (see buffers).
- Block non-specific binding sites with normal goat serum diluted 1/10 in TBS for 10 min at room temperature.
- Pour off the serum and remove excess serum tapping the slide edges on a paper towel.

Place 50-100  $\mu$ l of reconstituted anti-ISAV Mab onto the tissue sections (the volume added will depend on the size of sample to be covered) and incubate for 60 min at room temperature in a humid chamber.

Use appropriate controls i.e. known positive tissue as a positive control and uninfected tissue as a negative control; these should both be incubated with the reconstituted Mab and PBS separately.

Wash slides three times with TBS.

Add goat anti-mouse IgG biotin conjugate (1/50 in TBS) to the slides for 60 min

Wash slides three times with TBS

Add streptavidin-horseradish peroxidase (1/50 in TBS) to the slides for 60 min

Wash slides three times with TBS

To visualise the reaction, incubate the slides for 10 min with DAB solution (see buffers) or with a commercially available True Blue staining kit following the manufactures instructions.

If stained with DAB stop the reaction by immersing the slides in tap water and counter-stain them with haematoxylin for 3-4 min.

Rinse in tap water for 10 minutes.

Dehydrate the slides in 70% ethanol (3 min), 100% ethanol (5 min), xylene (2 x 5 min)

Mount the slides with Pertex and leave in fume cupboard to set.

Examine tissue under a light microscope –cells infected with the virus appear golden brown in colour when stained with DAB.

## Buffers

Phosphate buffered saline (PBS)

0.02M Phosphate, 0.15M NaCl      pH adjusted to 7.2 with HCl

$\text{NaH}_2\text{PO}_4 \cdot 2\text{H}_2\text{O}$                       0.876g/l

$\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$                       2.56g/l

NaCl                                        8.77g/l

Tris buffered saline (TBS)

Trisma base                              2.42g

NaCl                                        29.24g

Dissolve in approximately 900 ml distilled water, adjust pH to 7.2 using HCl and make up to 1 litre

10% (v/v) Hydrogen peroxide in methanol

Add 1ml  $\text{H}_2\text{O}_2$  (30% v/v solution) to 9 ml methanol

3,3'-Diaminobenzidinetetrahydrochloride (DAB)

Dissolve one 10mg tablet DAB in 6.67mls TBS

Place 0.5 ml aliquots of the solution into bijoux bottles, store at  $-20^\circ\text{C}$ .

For use add 5mls TBS and 0.1ml 1 %  $\text{H}_2\text{O}_2$  to 0.5 ml aliquot

NB. DAB is a possible carcinogen